


LE420

For demanding samples such as waste water, emulsions,
oily samples; refillable electrolyte

pH Ref °C ISM



Sensor type combined pH electrode

Order number 51340332

Specifications

pH range 0...14

zero point (pH) 7.00 ± 0.25

temperature range 0...80 °C

temperature probe -

shaft material glass

membrane glass HA

membrane resistance $< 600\text{M}\Omega$

reference system Ag / AgCl

reference diaphragm movable PTFE sleeve

reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head -

cable length 1.0 m

connector BNC

cell constant (cm^{-1}) -

cell type -

indicating metal -

Storage

3M KCl (51350072)