


LE422

For small samples in vials and test tubes

pH Ref °C ISM



Sensor type	combined pH electrode
--------------------	-----------------------

Order number	30089747
---------------------	----------

Specifications

pH range	0...14
zero point (pH)	7.00 ± 0.25

temperature range	0...80 °C
temperature probe	-

shaft material	glass
membrane glass	U
membrane resistance	< 250MΩ

reference system	Ag / AgCl
reference diaphragm	ceramic
reference electrolyte	gel
bridge electrolyte	-

electrode head	-
cable length	1.0 m
connector	BNC

cell constant (cm ⁻¹)	-
cell type	-

indicating metal	-
------------------	---

Storage

3M KCl (51350072)
