


LE438

For general applications in the laboratory
Robust and maintenance-free

pH Ref °C ISM



Sensor type	combined pH electrode
Order number	51340242
Specifications	
pH range	0...14
zero point (pH)	7.00 ± 0.25
temperature range	0...80 °C
temperature probe	NTC 30 kΩ
shaft material	POM
membrane glass	U
membrane resistance	< 250MΩ
reference system	Ag / AgCl
reference diaphragm	polyester
reference electrolyte	gel
bridge electrolyte	-
electrode head	-
cable length	1.0 m
connector	BNC
cell constant (cm ⁻¹)	-
cell type	-
indicating metal	-
Storage	
FRISCOLYT-B Electrolyte (51350076)	