

# GLPWAY PRINTER



## SPECIFICATION

Models	GLP Way				
Part NO	SGP - 011				
Security	3-levels of Password protection				
Printer unit	Dot matrix impact printer				
Ribbon cartridge	Exchangeable,Black ink				
Print format	5x7,24-character line length				
Paper roll	Standard paper 57.5mm(W),max.Ø 58mm				
Print quality	Thermally stable and light fast paper(GxP),(Non thermal paper)				
Print unit durability	1.5 million lines				
Power	230V AC,50/60 Hz				
Interface	RS232, 9Pin male connector				
Display		4.3"Colour LCD touch(resistive)			
Enclosure	ABS plastic				
Altitude	Up to 4000m				
Operating temperature	0-40°C Max,.80% at 31°C,linear decreasing relative humidity to 50% at 40°C				
Over voltage category	Class II				
Included in delivery	Power cable (country specific),interface cable(RS232 - RS232 only),paper roll(1pcs),Ribbon cartridge(1pcs),USB mouse(1pcs)				
Dimensions	185mm x 260mm x75mm(WxDxH)				
Net weight		600g (paper roll included)			



## Attributable

Create As many as 100 users, Let Every Operator log in with their Own User Name, Employee Code# & Password. The Documentation will bear the respective User Name.

09:09:14 01/11/2022

USERNAME **SCOUTER**

PASSWORD **\*\*\*\*\***

< >



## Legible

Eliminate Hand written entries and transcription errors. The Dot matrix print outs are legible & Clear.

SCOUTER PAPERWEIGHTS  
PRIVATE LIMITED  
Chennai

Ball, Make: Sartorius  
Ball, Mod: AR50245-CR  
Ser, No: 127990134  
S-W, Rev: 10P-23-07-00-WT  
Emp-ID: 100-ESP-481  
Batch, Number: 1900  
Lot No: 1101-0455  
Batch, ID: 011234  
AR NO: 101-0077  
Date: 11/07/2023  
St. Time: 11:51:00  
Wt: 10.00000g  
Wt: 10.00000g  
End Time: 11:51:07  
Op: Name: R. SURESH K. SURESH  
Sign: 1



## Original

Impact Printing will not fade away over time like Thermal Printing. So whatever the application & Documentation, All Original Documentation are safe for later reference / Verification.



## Accurate

The documentation completely reflects the true observation. This means an honest, accurate, & thorough representation of facts describing the conduct of the study.

Date & Time Log

No: 1

Note: Time drift Adjustment

Old date & time: 10/Jul/2023 12:33:13

New date & time: 10/Jul/2023 12:35:13

Cancel Print



## Security

Control unauthorized changes of settings Including Date & Time. It implies that standard privileged users can not access the Date & Time settings. Even if a privileged user accesses the menu and makes any changes to Date & time Settings, it will be appearing in the Audit Trail log file.



## Traceability

Add Alpha Numeric sample Identifiers to all your Samples as you Desire!.

09:12:43 GENERAL WEIGHTMENT COUNT 1

Register weight Batch No ba-123

WEIGHT Lot No lt-34

10.00000 g Batch Id ba-2

AR NO 01-987123456007

REGISTER NEXT WEIGHT

LINE FEED Tare Balance PRINT UNIT

Back MC