

## **SGP-04 Printer**



### **DOT MATRIX PRINTERS IDEAL FOR PLATFORM SCALES / LAB BALANCES / OTHER INSTRUMENTS**

- ⚙️ Eliminate Transcription errors Handwritten entries in Log books /Journal notebooks are time consuming and a source of transcription errors. Records on the normal printer paper helps to avoid such mistakes. The printouts can be archived and at any time provide proof of the original weighing data.
- ⚙️ Impact printing will not fade over time like thermal printing. So whatever the application, important results are safe for later reference.
- ⚙️ Durable & easy to use printer
- ⚙️ Every Individual Data received from the Platform scale/ Device will be printed.

Printer unit	Dot matrix impact printer
Ribbon cartridge	Exchangeable,Black ink
Print format	5x7,24-character line length
Paper roll	Standard paper 57.5mm(W),max.Ø 58mm
Print quality	Thermally stable and light fast paper(GxP),(Non thermal paper)
Print unit durability	1.5 million lines
Power	230V AC,50/60 Hz
Interface	RS232, 9Pin male connector